



king electronic tablet counting machine model TB2

for hospital dispensing



TB2
 Machine size: 39 x 24 x 28in.
 100 x 61 x 71cm.
 Voltage: 1 phase
 Operator: One
 Changeover time: 5 mins.
 Cleaning time: 5 mins.
 Power Consumption: 1/2 H.P.
 Net Weight: 200lbs. (90kg.)
 Gross Weight: 300lbs. (135kg.)
 Case Size: 39 x 24 x 27in.
 99 x 61 x 68cm.



Counts batches 1-1,000.
 It will handle any size or shape of tablet or capsule up to 16mm.
 Solid-state electronics.

Output is as follows:

Tablet	output per min.	Capsules per min.
5 mm.	190in.	No. 1 1000
6.5	250	2 1100
9	350	3 1200
10	400	4 1300
11.5	450	5 1400

TECHNICAL DATA TB2

The tablets are placed in a fibre glass hopper. A variable speed vibrator feeds the tablets over a dust-removing grid on to a rotating glass disc. The disc orientates the tablets and a single row drops off the edge of the disc through a light beam and photocell assembly. The pulses from the p/e cell are fed into an electronic batch counter. The count is selected on rotating switches. Any batch can be selected from 1-1000 (3-decade) or 1-10000 (4-decade). The tablet feed is continuous and at the end of the batch the tablets are directed into alternate containers.



Handing ApS
 Nordkranvej 5-9, Vassingerød
 DK - 3540 Lyngø
 Tel.: +45 48160166
 www.handing.dk

King TB2

Dette dokument er vejledende.
 Der tages forbehold for fejl.
 This document is only guiding.
 We make reservation for possible errors.

Tegn. Godk.TN

Lager nr.

F423